

# ToxReader



## Interpret ToxWipe Results with Confidence

Compatible with ToxWipe 6, 6+ and 6 Elite oral drug testing devices. Scan a unique device QR code before each test. Stores and prints results for easy access records. ISO 13485:2016 and CE certifications.

CPU	NXP4330Q/Cortex-A9 Quad Core /1.4GHz
Display	4.3" Touch Screen(480(RGB)* 800 Pixels)
Operating System	Android 5.1.1
Dimensions	218*90*56.6mm
Weight	442g
Battery	7.4V/2100mA Lithium Battery
Battery Life	350 tests (with printing) 1000 tests (without printing)
Camera	5 Mega Pixel
RAM	1GB DDR3, 16GB eMMC
Storage	100,000 tests records
GPS	GPS C/A Code
Printer	Thermal Printer/External Printer
Internet	WIFI and 3G

**1-800-559-2490**  
**www.oranoxis.com**

